L.N.S.[®] GT

L.N.S.[®] GT - STANDARD Y LARGO
L.N.S.[®] GT - STANDARD AND LONG

surgival
GRUPO COSÍAS

Índice / Index

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[Schipper et al. J Bone Joint Surg (Br) 2004;86-B;86-94].

Fractures of the upper third the femur are currently treated using different models of intramedullary trochanteric nails. Trochanteric intramedullary nailing is a demanding technique, sometimes entailing unexpected problems.

A satisfactory surgical result involves various factors and depends on:

- Good preoperative and surgical planning.
- Acceptable fracture reduction.
- Good sagittal and axial viewing of the fracture.
- Anatomical and spatial knowledge of the relationship between the greater trochanter and the femoral diaphysis.
- Reliable instruments.
- Strictly followed surgical stages.
- Post-operative care.

The most common post-operative complication is cut-out of screws from the femoral head.

“About 80% of the operations in patients in whom cut-out of the implant occurred were reported as being moderate or very difficult, the reduction non anatomical but acceptable and the screws in suboptimal position”.

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[Anteroposterior x-ray study of the proximal third of both femurs.](#)

- [Study of the contralateral femur with the leg in neutral rotation may give valuable information about the normal cervical-diaphyseal angle.](#)
- [The type of implant to be used depends on the type of fracture, which is best seen by applying distal traction to the fractured leg and keeping it in neutral rotation while the x-ray study is performed.](#)

Before starting the operation, be sure the following are at hand:

- FOR THE STANDARD NAIL:
 - Flexible reamers.
 - A SHS (a very narrow femoral canal or an irreducible fracture might require open reduction and a sliding hip screw implant).
 - A L.N.S.[®] GT-LONG implant (a preoperatively unseen subtrochanteric crackline or fracture or an iatrogenic fracture might need a long implant locked in the distal femur).
- FOR THE LONG NAIL:
 - Flexible reamers.

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Place the patient on the traction table on his back, keeping his pelvis parallel to the floor.

Adduce the fractured leg slight, with the hip flexed 10° to 20° and the trunk bent towards the opposite side, maintaining this position to allow easier surgical access to the greater trochanter.


Fig. 1

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The unaffected leg should be abducted as far as possible, or abducted with the hip flexed (Fig. 1).


Fig. 2

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Place an image intensifier in such a way as to provide clear images in the anteroposterior and axial plane (Fig. 2).


Fig. 3

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The angle of inclination of the image intensifier in the axial view which shows the femoral shaft, neck and head aligned in a straight line, matches the anteversion angle of the femoral neck (Fig. 3-4).

Fig. 4

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Fracture reduction should be as anatomical as possible. An implant should never be placed on an unacceptably reduced fracture.


Most pertrochanteric fractures are reduced by longitudinal traction and varying amounts of internal rotation of the lower limb. Reduction must be correct in both planes.

Very unstable fractures tend to have an angulation seen in lateral view as posterior sag (Fig. 5). Manual postero-anterior pressure on the apex of the angulation (Fig. 6) during surgery corrects this deformity (Fig. 7).

Sometimes, both bone fragments are parallel, but the distal fragment is posterior (Fig. 8). If manual postero-anterior pressure (Fig. 6) reduces the fracture (Fig. 7), the surgical field is ready and the operation can begin. Insert the assembled L.N.S.®GT in the femoral canal.

Pushing the whole construct upwards (Fig. 9) reduces the fracture and allows insertion of the neck screws in the proper position*.

If neither of these operations produces acceptable reduction, open reduction will be required.

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The surgical field should extend from the base of the thorax to the knee joint, to allow free access to the greater trochanter and so that, if necessary, a long nail can be locked in the distal end of the femur without the need for preparation of a larger surgical field.

(*) Intraoperative closed reduction and percutaneous fixation of displaced femoral neck fractures with the trochanteric gamma nail.
Christos Garnavos. Tips & Techniques. Orthopedics. (March 2003)

TÉCNICA QUIRÚRGICA CLAVO L.N.S.[®] GT - STANDARD SURGICAL TECHNIQUE L.N.S.[®] GT - STANDARD

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To facilitate instrumentation of the nail, it is important that the cutaneous incision be placed 5 cm proximally to the profile of the greater trochanter. After fracture reduction, with the help of the C-arm and a K-wire, two skin reference lines are drawn:

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
- b) the lateral projection of the tip of the greater trochanter as seen in the A-P view (Fig. 11a y 11b). The cutting point of these two lines, locates the tip of the trochanter in both planes.

Fig. 11a

Fig. 11b

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Start the incision 4 to 5 cm proximally to the apex of the trochanter. Extend the incision cranially for 2 to 5 cm in line with the lateral axis of the femur (Fig. 12). When inserting the L.N.S.[®]GT - STANDARD nail, the skin incision can be made slightly anterior to the lateral femoral axis, in order to give better access to the anterior third of the greater trochanter.

Cut down through the subcutaneous tissues. Opening the gluteus maximus fascia in line with its fibers, allows you to reach the upper edge of the greater trochanter with the tip of your index finger (Fig. 13).

It is easy to identify the postero-superior corner of the greater trochanter with your fingertip (Fig. 14).


Fig. 12

Fig. 13


Fig. 14

2. Portal de entrada femoral (trocantereo y diafisario) Femoral entry point (trochanteric and diaphyseal)

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
The entry point is located at the junction of the anterior third and posterior two thirds as seen in the AP view, at or slightly lateral to the upper edge of the greater trochanter (Fig. 15). In the lateral view the entry point is centered along the axis of the femoral canal, or slightly anterior to it (Fig. 16).

Fig. 15

Fig. 16

APERTURA FEMORAL (DOS POSIBILIDADES TÉCNICAS) OPENING THE FEMUR (TWO POSSIBLE TECHNIQUES)

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a) Percutaneous technique

Insert a 2.5 mm guide-wire and check the right entry point over the greater trochanter with the C-arm. Slide the cannulated curved awl along the guide wire and with manual rotary movements create the entry hole (Fig. 17). Remove the curved awl keeping the guide-wire in the femoral canal. Slide the soft tissue protection sleeve along the guide-wire. Remove the inner sleeve and slide along the guide-wire the manual single-pass straight reamer (Fig. 18) or a power-driven 16-17 mm flexible reamer to two centimeters below the level of the lesser trochanter.

With osteoporotic bone or comminuted greater trochanter, the femur can be opened easily by inserting the manual single-pass straight reamer through the guide-wire (Fig. 18).

b) Standard technique (incision 4-7 cm)

Depending on the surgeon's criteria, in especially difficult cases, opening of the greater trochanter can be done with a standard awl.

Fig. 17


Fig. 18

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Attach the nail to the nail-guide (Fig. 19). Insert the locking bolt into the nail-guide and use the T-handle locking-wrench to secure the bolt to the nail (Fig. 20), and tighten well. Screw the small threaded driver to help insertion of the nail, if necessary (Fig. 21).

Fig. 19

Fig. 20

Fig. 21

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The L.N.S.® GT - STANDARD nail is introduced manually mounted in its guide with rotary movements and x-ray control. An extension can be added to the guide for easier insertion or removal. (Fig. 21).

When the nail to be implanted has been selected according to the pre-operative study, mount it in the guide and check the guide-hole settings of the sleeves and drills for locking (on the instrument table) (Fig. 22 and 23).

Fig. 22


Fig. 23

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When the nail has been fully introduced, make slight rotary movements to see the notches in the hole for the dynamic neck screw. This maneuver shows you the future situation of the neck-screw, by placing a K-wire over the anterior surface of the hip and centering it over the round notches of the nail.

Fig. 23b

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
The nail should previously have been arranged in anteversion to the hip. An instrument is available which allows the screws to be positioned using the Kirschner wire. This position should coincide with the direction of the neck to ensure that the nail is placed in the desired anteversion (Fig. 24).

Fig. 24

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The length of nail's insertion can be marked with a guide-wire passed through a hole in the nail-guide and it is anchored into the lateral cortex or into the soft tissues. The end of the K-wire does shows the nail-guide-nail union (Fig. 23b).

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
The system has a Ø 3.2 mm Kirschner wire with a drill point and a threaded area allowing greater cutting efficiency and reduction of the wire bending to a minimum. The wire is also identified by the color yellow, like all the instruments related with the compression screw.

The 3.2 mm guide-wire can slide around the smooth, round femoral surface and enter eccentrically through the nail hole. To avoid this complication, the straight awl and its sleeve are inserted through the compression-screw reamer sleeve. An incision is made at the skin contact point allowing the construct to be pushed through the muscle. The awl and sleeves are pushed through the Vasturs lateralis until they reach the lateral femoral cortex. Turn the awl manually to make a small notch in the bone as a starting point for the lag-screw 3.2 mm guide-wire (Fig. 25).


Fig. 25

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Remove the straight awl and its sleeve and replace it with the internal sleeve with the 3.2 mm K-wire.

It is very important to keep the compression screw sleeve tight up to the femoral cortex when inserting the 3.2 mm K-wire in the femur.

The tip of the K-wire is a drill point to facilitate intraosseous insertion without deviation or deflection (Fig. 26).


Fig. 26
Punta broca roscada
de la aguja.
Threaded drill tip
of the wire

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
The 3.2 mm guide-wire should be moved under C-arm control, checking the guide-wire position in both planes. The sharp point of the 3.2 mm guide-wire should reach or be within 2-3 mm of the subchondral bone (Fig. 27).

Fig. 27

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In anteroposterior view, the needle should be in the bottom half of the femoral neck. In axial view, it should be centered or slightly to the back of the femoral head and neck. The guide wire should not be in the top of the femoral head and neck, because that would increase the possibility of cut-out.

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Remove the guide wire sleeve and introduce the length meter (Fig. 28), which should be tightly inserted (Fig. 29). It will indicate the screw's length (Fig. 30). If a proximal screw is used, its guide wire must be placed after placement of the neck screw.

Fig. 28

Fig. 29

Fig. 30

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The length of the compression screw is marked with the movable stop to the extent corresponding to the drill's graduated scale (Fig. 31). It is recommended that the stop be set 5.0 mm below the measurement obtained with the length meter, if the guide wire reached the subchondral bone (measured length: 100 mm; drill set to 95 mm).


Fig. 31

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Under C-arm control, the drill is passed along the guide wire and inside the femoral head until the drill's movable stop touches the sleeve (Fig. 32).

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The drill point of the dynamic screw and the screw may jam or touch the metal edge of the nail hole. In this case, the guide wire must be removed and the drill point and screw be reinserted without the wire. .

Fig. 32

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The compression screw inserter has a bar inside which threads to the screw to hold it (Fig. 33). The assembly is passed through the sleeve and moved to the desired position in the femoral head (with C-arm control) (Fig. 34).

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If smooth, moderate interfragmentary compression is wanted, turn the nut clockwise. Compression begins when the wheel makes contact with the heuristic sleeve (Fig. 34).


Fig. 33


Fig. 34

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The system has a Ø 2.5 mm Kirschner wire with a drill point and a threaded area allowing greater cutting efficiency and reduction of the wire bending to a minimum. The wire is also identified by the color black, like all the instruments related with the anti-rotational screw.

Place the sleeve of the anti-rotation screw drill tip with the internal sleeve for the 2.5 mm guide wire through the anti-rotation hole of the nail guide (Fig. 35). Make the incision at the point of skin contact and push the sleeves until they make contact with the femoral cortex.

The sleeves and awl are pushed through the Vastus lateralis until they make contact with the lateral femoral cortex. Turn the awl until it makes a small hole in the femoral cortex to serve as the drilling start point for the 2.5 mm wire.

Remove the straight awl and its sleeve and replace it with the internal sleeve with the 2.5 mm guide wire.

It is very important to keep the compression screw sleeve tightly pressed on the femoral cortex when beginning insertion of the 2.5 mm guide wire in the femur.

The tip of the guide wire is designed as a drill point to facilitate intraosseous approach without deviation or deflection.

With x-ray control, introduce the guide wire to the desired length in the femoral head.

Fig. 35

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Remove the guide wire sleeve and insert the length meter to obtain the screw length. It is recommended that the screw should be about 20 mm shorter than the neck screw (Fig. 36).

Fig. 36

6.3. Perforación femoral para el tornillo proximal / Proximal screw femoral drilling

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Follow the same steps as in Section 5.3. (Fig. 37 and 38).

Fig. 38


Fig. 37

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Remove the length meter and pass the 6.0 mm drill tip with the guide wire to the selected length. Insert the screw through the sleeve with the hexagonal screwdriver (Fig. 39). Passing the drill tip and placing the screw should be done with C-arm control.

Fig. 39

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Select the guide hole for making static or dynamic distal locking. Through this, pass the 4.5-mm or 3.5-mm screw drill tip, for 11-mm or 10-mm diameter nails, respectively.

At the cutaneous contact point, make the incision allowing the sleeve to be pushed through to the femoral diaphysis. Insert the calibrated drill bit through the sleeve and perforate both cortices. Remove the drill bit and place the length meter to obtain the screw length.

When you have selected the screw, put it in the hexagonal screwdriver (Fig. 40) and place it through the sleeve.

This whole process should be C-arm controlled.

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After placing the femoral head compression screw, without removing the screw inserter, adduce the hip to make insertion of the end-cup easier. Use the T-handle locking-wrench to remove the locking bolt carefully, holding the nail guide still manually so that the joint is not broken: keeping the dynamic screw inserter in place makes keeping the connection easier.

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The end-cap is fixed to its inserter shaft with the locking bolt (Fig. 41) and the T-handle locking-wrench. The assembly is inserted through the top hole of the nail guide until it makes contact with the proximal end of the nail (Fig. 42) then screwed on.

Fig. 41

Fig. 42

TÉCNICA QUIRÚRGICA L.N.S.[®]GT - LARGO **SURGICAL TECHNIQUE L.N.S.[®]GT - LONG**

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Steps 1, 2, 3, 5, 6 and 8 of the L.N.S.[®]GT - STANDARD nail technique should be repeated in the same way for the L.N.S.[®]GT - LONG nail system. These are the specific steps:

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The L.N.S.[®]GT - LONG nail allows the blunt-tipped intramedullar mill guide to be passed through the nail's 4.5 mm cannulation. So there is no need to change the blunt-edged guide for a sharp-ended one (Fig 43). The diaphyseal femoral canal is milled with the flexible drill.

Fig. 43

2. 2. Inserción del clavo L.N.S.[®]GT - LARGO / L.N.S.[®]GT - LONG nail insertion

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The L.N.S.[®]GT - LONG nail has an anteroposterior radius of curvature to fit anatomically most femur deflections. The 13° anteversion of the holes for the neck screws is preset in the angle of anchorage of the nail with its guide (Fig. 44).

Fig. 44

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To facilitate insertion of the L.N.S.® GT - LONG nail, the guide with the nail in place is first turned 90° so that the anteroposterior deflection of the nail is seen entirely in the medial-lateral plane (Fig. 45).

Insertion is begun maintaining this position and the guide is pushed distally with small rotational movements, until the point of the nail reaches the femoral isthmus.

To drive the nail further after this point, the nail guide is turned progressively in the backwards direction until it reaches the plane parallel to the floor (Fig. 46 and 47) and pushed distally.

Sometimes it is necessary to strike the guide with a hammer to insert the nail. In these cases, it is necessary to be careful and check the progress of the nail after every hammer blow. If the hammer blows do not push the nail in, it is recommended that the nail be removed and the medullar canal be drilled with a drill of greater diameter.

Fig. 45

Fig. 46


Fig. 47

3. Enclavado distal / Distal locking “Técnica de manos libres” / “Free hand technique”

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Distal locking can be performed with one or two locking screws. The L.N.S.[®] GT - LONG nail offers these options:

1. Single static locking screw through the upper distal round hole.
2. Single static locking screw through the proximal zone of the lower oval distal hole.
3. Single dynamic locking screw through the distal zone of the lower oval distal hole.
4. Double static locking screws through positions 1 and 2.
5. Double locking screws through positions 1 and 3.

Single distal locking screws are mainly indicated for single or segmental fractures located in the upper half of the femur.

- AO/OTA 31A-3. Transtrochanteric and obliquity reverse fractures.
- High and Low Subtrochanteric Fractures.
- AO/OTA 31A-2 and 31A-3 fractures combined with non-comminuted fractures of the upper middle half of the femur.
- Impending pathologic bone lesions and pathologic fractures of the upper third of the femur.

A single dynamic locking screw is indicated in cases where it is desirable to allow compression forces acting between the bone fragments (AO/OTA 31A-3 and Subtrochanteric fractures).

AO/OTA 31A-2 unstable fractures can be treated with the GT - LONG nail, with no distal locking screw.

Double distal locking is mainly indicated in segmental comminuted fractures, or when the distal fracture is located in the lower half of the femoral shaft. Due to the wider femoral canal, the distal fragment is best stabilized by using two distal locking screws:

- AO/OTA 31A-2 and 31A-3 fractures combined with comminuted upper or lower middle half femoral shaft fractures.
- AO/OTA 31A-2 and 31A-3 fractures combined with non-comminuted lower middle half femoral shaft fractures.

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The surgical technique to place the femoral neck-head screws or the compression-screw in the optimal position into the femoral head is the same as in the L.N.S.[®] GT - STANDARD nail for AO/OTA 31A-1, 31A-2 and 31A-3 fractures.

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Recommendations:

- Rehabilitation should be commenced early.
- Active and passive movements of the knee and ankle should be carried out in the first 24 hours after the operation.
- Full weight bearing with external help should be begun as soon as possible.

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Through small cutaneous incisions, remove the screws of the distal locking with the hexagonal screwdriver.

Placing a 2-mm guide wire through the screws of the neck may assist in their extraction with the dynamic screw inserter and hexagonal screwdriver.

To facilitate extraction of the end-cup and nail, the hip may be adducted. The cup is removed with the T-handle locking-wrench connected to its insertion shaft (Fig 47). Screw the locking bolt on to the nail and use a slap hammer to remove the nail from inside the femoral canal.

Fig. 47

BIBLIOGRAFÍA RECOMENDADA / RECOMMENDED BIBLIOGRAPHY

1. Failure of femoral head fixation: a cadaveric analysis of lag screw cut-out with the gamma locking nail and AO dynamic hip screw.
R.C. Haynes, et al. Injury(28) No5,6; 337-341:1997.
2. Use of an Intramedullary Hip-Screw Compared with a Compression Hip-Screw with a Plate for Intertrochanteric Femoral Fractures. A prospective, randomized Study of One Hundred Patients.
Hardy, Dominique et al. J Bone Joint Surg. (Am) 1998 ;80-A ; 618-30
3. Slotted Intramedullary Hip Screw Nails Reduce Proximal Mechanical Unloading.
Dominique C.R. Hardy, and Konstantin Drossos. Clin Orthop 406, 176-184:2003
4. A Biomechanical Evaluation of the Gamma Nail.
Shepherd F. Rosenblum, et al. J Bone Joint Surg. (Br) 1992; 74-B:352-7.
5. Biomechanical Analysis of Location of Lag Screw of a Dynamic Hip Screw in Treatment of Unstable Intertrochanteric Fracture.
Chi-Chuan Wu et al. J Trauma 41(4);699-702:1996.
6. Intramedullary nails: some design features of the distal end.
Wang CJ, Brown CJ, Yettram AL, Procter P. Med Eng Phys. 2003 Nov;25(9):789-94.
7. Intramedullary femoral nails: one or two lag screws?. A preliminary study.
Wang CJ, Brown CJ, Yettram AL, Procter P. Med Eng Phys. 2000 Nov;22(9):613-24.
8. Intramedullary nails with two lag screws.
Brown CJ, Wang CJ, Yettram AL, Procter P. Clin Biomech (Bristol, Avon). 2004 Jun;19(5):519-25.
9. Percutaneous Treatment of Peritrochanteric Fractures Using the Gamma Nail.
Carlo Bellabarba, Dolfi Herscovici, Jr and William M. Ricci. Clin Orthop 375;30-42;2000.
10. Locking and Compression of the Lag Screw in Trochanteric Fractures is not Beneficial. A prospective randomized study of 153 cases.
Sernbo I, Johnell O, Gardsell A. Acta Orthop Scand, 1994 Feb;65(1):24-6
11. Geometric Mismatch of the Gamma Nail to the Chinese Femur.
K.S. Leung; P. Procter, B. Robionek, and K. Behrens. Clin Orthop 323:42-48;1996.
12. Multicenter trial of modified gamma nail in East Asia.
K. S. Leung et al.
Possibilities for Improving Fatigue Properties of Interlocking Screws of Solid Tibial Nails. A Mathematical Model with Practical Conclusions.
Clin Orthop 323;146-154;1996.
13. Thread Types of Locking Screws Concentration Factor and Edge Fiber Stress of Different Thread Types in Locking Screws. A Mathematical Approach Der Unfallchirurg 2002; 105: 140-146
Holper B, Tschegg EK, Stanzl-Tschegg S, Gabbler C. Unfallchirurg 2002 Feb;105(2):104-6.
14. Fatigue Stability of Locking Screws for Unreamed Tibial intramedullary Nailing.
Gabbler C, Tschegg EK, Greitbauer M, Stanz-Tschegg S, Wozasek GE, Laube W. Vecsei V. Unfallchirurg 1999 Jan;102(1):29-34.
15. Fatigue Strength of Locking Screws and Prototypes Used in Small-Diameter Tibial Nails: A Biomechanical Study.
Gaebler C, Stanz-Tschegg S, Heinze G, Holper B, Milne T, Berger G, Vecsei V. J Trauma 1999 Aug;47(2):379-84.
16. Bending Strength and Holding Power of a Prototype Tibial Locking Screw.
Lin J, Hou SM. Clin Orthop 2002 Oct;403:232-9
17. Bending Strength and Holding Power of Tibial Locking Screws.
Lin J, Lin SJ, Chiang H, Hou SM. Clin Orthop 2001 Apr;(385):199-206.
18. Mechanical Strength, Fatigue Life, and Failure Analysis of Two Prototypes and Five Conventional Tibial Locking Screws.
Sheng-Mou Hou, Jaw-Lin Wang, and Jinn Lin. J Orthop Trauma. 2002; 16 (10). 701-08.

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

9 Aguja Kirschner punta broca Ø 2,5 mm
Ø 2.5 mm K-wire with drill bit tip
Ref. C6235980 3 unidades / 3 units

Set 2 instrumental L.N.S.® GT L.N.S.® GT instrument set 2

Ref. C6236610

1 Punzón iniciador L.N.S.® GT
L.N.S.® GT starter awl
Ref. C6236040


2 Punzón recto
Straight awl
Ref. C6236085

3 Fresa iniciadora Ø 16 mm
Ø 16 mm starter drill
Ref. C6236045

10 Aguja Kirschner punta trocar Ø 2,5 mm
Ø 2.5 mm K-wire with trocar tip
Ref. C6236025 3 unidades / 3 units

11 Maza extractora
Extractor mallet
Ref. C6236105

12 Vaina-guía para fresa iniciadora
Guide-sleeve for starter drill
Ref. C6236050

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