Técnica Quirúrgica
Surgical Technique

GENUTECH®

Sistema total de rodilla primaria Genutech®
Genutech® total primary knee system

INSTRUMENTACIÓN PARA CORTE FEMORAL REDUCIDO
INSTRUMENTATION FOR REDUCED FEMORAL CUTTING

 **surgival**
GRUPO COSÍAS

Índice / Index

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A front central incision is used, with internal parapatellar arthrotomy, and external dislocation of the patella which, in the proximal direction, extends up to the middle third of the quadriceps tendon, and distally, to the tibial tuberosity.

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
After the approach, the following order is advised:

1. Tibial cutting
2. Front and back cutting of the femur
3. Distal cutting
4. Measurement of 90° flexion and extension spaces
5. Complete femoral cutting (bevels)
6. PS/NPS selection
7. Intercondylar box cutting (in case of PS)
8. Placement of femoral testing prosthesis
9. Preparation of tibial keel
10. Action on patella
11. Verification of mobility and stability with testing implants
12. Placement of final prosthesis

Closure in layers, optionally leaving drainage points, at the discretion of the surgeon.

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
The design allows for osteotomy with the extramedullary or intramedullary guide, or both being used together.

The tibial cutting template (left or right) is attached to the extramedullary rod which, in turn, distally secures adjustable clamps to the ankle.


If the extramedullary guide is going to be used, in the proximal area, a short auxiliary rod is assembled, equipped with two spicules attached to the area of the tibial spines once the position of the extramedullary rod has been decided.

If the intramedullary guide is going to be used, the tibial intramedullary canal is opened previously with the awl and the drill bit (the entry point is located between the tibial spines). Through the opening we have made, we will insert the intramedullary shaft, assembling the extramedullary rod in the proximal area.

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In cases where the tibia is pointing inward, it is advised to make the entry point a little more external, and if it is pointing outward, a little more medial. The intramedullary guide is placed on the upper tray of reduced cut femoral set 1.

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The probe is installed in the cutting guide's housing, and with the “10” short probe end, we will seek support on the least affected plateau, temporarily attaching the cutting guide at this level. After this, with the “0” long probe end, we will find the most depressed area of the tibia on the most affected plateau, the following two situations then being possible:


- **That the “0” end does not make contact**, which means we must lower the tibial cutting guide until making contact. At this cutting guide level, we ensure that we make a 1 mm resection below the most worn area.
- **That the “0” end does not fit**. In this case we will block the cutting guide.

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The cutting guide is attached to the tibia with four headless pins, inserted through the most distal orifices (the line marked with the "0").
Attachment is completed with the option of one or two pins with oblique-shaped and converging heads.

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
With the cutting guide firmly attached to the tibia, the intramedullary and/or extramedullary guide system is removed, and osteotomy of the tibia's epiphyseal area is performed.

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
Once the cut is made, to extract the guide, the oblique, headed pin is extracted. Once the tibia fragment is removed, any peripheral osteophytes are cut with a gouge, and the tibial tray is measured with the templates.


After this point, femoral instrumentation is begun.

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
It is recommended not to remove the four headless pins until the distal and posterior femoral cuts have been carried out, since they may be useful if a recut must be done (see Annex 2).

Set femoral 1 corte reducido. Bandeja superior
Reduced cutting femoral set 1. Upper tray

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
With the aid of the "T"-shaped quick coupling handle, the intramedullary shaft is inserted through the perforation previously made, leaving enough outside the perforation to allow for installing the femoral sizing guide. Once the intramedullary shaft is inserted, the "T" shaped handle is removed.

Set femoral 1 corte reducido. Bandeja superior
Reduced cutting femoral set 1. Upper tray

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

With the knee at 90°, the size of the femur is measured by inserting the intramedullary shaft through the cannula of the femoral sizing guide, moving it along until making contact with the femur's distal condyle. This position, when compared to that of the intramedullary shaft, makes it possible to take a reading of the femur's size without errors due to the unsuitable placement of the sizing guide.

The depth of the probe stylus is adjusted, ensuring that the reading recorded on it coincides with the size shown on the sizing guide's distal surface. These readings are taken from the position occupied by the surgeon.

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
On the front and back cutting mask ("A"), with the size determined by the femoral sizing guide, the upper lateral-medial rail is mounted, with the probe stylus of the corresponding size, as well as the Varus-Valgus Angle Adaptor (3°, 5° or 7°), determined in the pre-operation telemetry, which is best suited to the patient's morphometry. Said adapter is mounted by sliding down the central slit of mask "A".

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
The most appropriate external rotation system is pre-mounted (Neutral 0°, 3° or 6°), taking into account if the knee to be operated on is the right or the left one. To this end, the "Y" turn system must be mounted on the tibial base plate with the attachment turret in the lower position, with the corresponding marking to the Left ("L") or Right ("R"), as seen from the surgeon's point of view.

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Turning the handle, the external "Y" rotation system is attached, mounting it in its highest position so as to make it easier later on to adjust the height of the "A" cutting mask on the femoral component.

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
Once the rotation system's position is set, its lugs are introduced into the lower orifices for "A" mask pins, creating final insertion pressure, which will ensure the attachment of the mask.

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
Once distal contact is made between, at least, the most prominent condyle and the "A" cutting mask, the external rotation system is unlocked by turning the handle to the left, and the system is compressed, ensuring subsequent condyle contact with the skids.

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To ensure that the system remains compressed, the handle is turned to the right to tighten and secure its position.

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
Remove the tibia pins (left as reference points in case a hypothetical recut is necessary) if, upon lowering the system, it interferes with these pins.

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Then, the mask is attached to the bone with pins after drilling a Ø3.2 mm hole.

It is recommended that the distal pin first be attached through the front orifice resting on the most prominent condyle, and then attach the oblique pins located on both sides of the mask.

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
Alternatively, one may also opt to properly place the “A” cutting mask to carry out front and back cutting by following the steps in Annex I.

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With the "A" mask attached, the rotation system is taken out, along with the lateral-medial rail, pulling vertically on the probe stylus' threads. We insert the intramedullary shaft a little further, if necessary, to prevent it from blocking the saw, so that it is possible to perform the front and back cutting comfortably.

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
With the aid of the previously made cut, the "B" distal cutting mask is attached to the "A" mask, securing its position with pins.

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Once the “B” mask is attached, the intramedullary shaft is removed, as well as the “A” mask, firstly extracting the distal pins to avoid collisions with the cutting saw.

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
When removing the intramedullary shaft, the “Angle Adaptor” located in the central slit runs the risk of falling, so it is recommended to hold it in place during this maneuver and also remove it later.

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
With the “B” Mask attached, the distal cut is made.

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Once the cut is made, the “B” mask is removed.

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

Having performed the tibial cut, and with the distal and posterior cuts made, the spacers are inserted in flexion and extension to evaluate limb alignment and ligament tension.

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If excessive ligament tension is observed, see Annex 2 to make “recuts”.

Set femoral 2 corte reducido
Reduced cutting femoral set 2

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On the Front Cut, and ensuring proper distal support, we attach the “Bevel Cutting Mask” corresponding to the size of the femur, securing it to the bone with the lateral pins.

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Do not use frontal pin positions, since they go through the bevels and would impede the cutting saw's passage (See the following Note).

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Once the masks are mounted, the Back Bevel Cutting is performed (recommended to maintain maximum distal support area), and finally the Front Bevel Cutting.

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

If you observe that the mask is not securely attached, it may be improved as follows:

- **Inserting frontal pins.** One must take the precaution of always inserting them in the part opposite to where we are making the bevel cuts in the condyle on which we are working.
- **By again mounting the “B” Distal Mask**, and attaching it with pins so that it also provides stability to the “A” Mask during cutting.

Set fémoro-patelar. Bandeja inferior
Patella-femoral set. Lower tray

Set PS. Bandeja superior e inferior
PS set. Upper and lower tray

2.10 Implantación del componente femoral de prueba / Implantation of femoral testing component

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
The masks are removed and the “NPS” femoral testing component is inserted, attaching it to the femur with two front-distal pins.

- If one chooses to implant a posterior stabiliser (PS), the NPS component is coupled with the saw guide to recut the intercondyle area. Later on, the PS femoral test component is attached.

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

With the femoral test component installed (either PS or NPS), the tibial tray is inserted, with the test insert of a pre-determined size and thickness, supported on the tibia's cutting surface.

With the handle which attaches the tray, the rotation of the tibia component is adjusted, verifying proper flexo-extension behaviour in relation to the femoral test component previously installed. The tray must be pointing approximately toward the ATT.

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
When it is clear that it is in a suitable rotation position, a mark is made between the tray and the bone which will act as a reference point, so that later on, once the test polyethylene is removed, the tray can be attached with two headed pins.

The keel-chisel of the chosen size is inserted as far as possible into the tibial tray. The keel handle is removed, and the test insert is attached to perform a final verification of the corresponding tibial and femoral test components.

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

After performing said verification, the test insert and keel are removed, deepening the central cylinder as is appropriate to house the tibial tray cover, according to whether it will be “standard” or “long”.

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With a prosthetic replacement for the patella, there are two aspects to take into account:

- Knowing the original height of the patella, it is recommended leave a minimum bone thickness to avoid fracture risk (12-14 mm).
- Avoid excess pressure which would result in a prosthetic patella being higher than the original. To this end, when the size of the prosthetic patella is selected, one must not only be guided by the diameter which best fits the bone resection surface, but also by the height of each patella component measurement:

Diámetro / Diameter	32 mm	34 mm	36 mm	38 mm	40 mm
Altura / Height	7 mm	8 mm	9 mm	10 mm	11 mm

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
The osteophytes and the synovial tissue surrounding the patella border are then removed, until the quadriceps and patella tendon are visible.

With the patella resection clamps, closed over the patella so that its frontal cortex rests on the bolts, and the turret is adjusted at 12 – 14 mm, the patella cut is made.

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The clamp's mobile arm is positioned so that the perforation of the prosthetic pivot is slightly medialised in relation to the centre of the patella. The perforation is performed with the specific drill bit.

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Once the test patella is in place, its movement is evaluated through the femoral intercondylar channel with flexo-extension movements.

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If, with the test (or prosthetic) patella, tension in the extensor apparatus increases, and with it, possibilities of subluxation and wear, a smaller patella component would be chosen to decrease its height.

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After having performed appropriate verifications with femoral, tibial and patellar test components, they will be removed and replaced with the final prostheses.


Once the prosthetic solution is installed, we will again verify that the proper function previously observed with the test components has been maintained, after which we will finally attach the tibial insert, screwing the front bolt to the tibial tray.

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
With the probe stylus supported on the front part of the femur, with the knee at 90° flexion, and with the entire assembly only introduced on the intra-medullary shaft, the alignment shaft must be aligned with the tibia load axis to obtain the mask's external rotation.

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Maintaining this position, it is recommended to first place the distal pin through the front orifice resting on the most prominent condyle, and then secure the oblique pins located on both sides of the mask. Once the mask is attached, we will make the front and back cuts.

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
If the spacers are introduced in flexion and extension, and it is judged that ligament tension is excessive, it may be due to the following causes:

- Does not enter in flexion or in extension.

It is recommended to make an additional 2 mm recut on the tibia (this is why it is not recommended to remove the tibial cutting guide attachment pins until reaching this point, since we will only have to insert the mask through the orifice line marked with “+2”). With it, we will lower the cutting guide's position by 2 mm from the “0” position upon which we initially made the tibia cut.

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Before performing the cut with the saw, the oblique headed pin(s) are again attached to avoid the movement that the saw's vibration may cause to the cutting guide:

- Enters in extension, but not in flexion.

It is recommended that the possibility of detaching the posterior femur capsule be evaluated, or to choose a smaller femoral size (especially in cases where the objective is to preserve the posterior cruciate ligament, since it cannot be detached).

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- Enters in flexion, but not in extension.

It is recommended that a distal femoral recut be made, of the thickness judged appropriate (taking into account the modification of the inter-line).

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The femoral "Distal Recut" is done by supporting the "Recutting Mask" on the front and distal levels.

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Once the distal recut is made, we return to the bevel position. Once we have reached this point, two situations may arise:

1ª. We have verified with the spacers (page 16) and, consequently, we do not have the bevels cut.

In this case, after the distal recut, we could use the bevel cutting masks to cut the edges and continue with the surgical technique.

2ª. That we have made the bevel cuttings, so we will again mount the "A" Cutting Mask on the "B" Distal Mask, again supporting it on the cuts present, and attaching it with oblique pins (located on the sides of the mask).

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Do not use frontal pin positions since they go through the bevels and would impede the cutting saw's passage.

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Once the entire assembly is attached the bevel cuts are made.

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If you observe that the mask is not securely attached, it may be improved as follows:


- Inserting frontal pins. One must take the precaution of always inserting them in the part opposite to where we are making the bevel cuts.
- Put pins on the "B" Distal Mask as well so that it provides stability to the "A" Mask while re-cutting bevels, in exceptional cases of extreme attachment deficiency.

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Clipping is done manually, first supporting the back part of the insert on the tibial tray and pressing both pieces on the front until clipping is achieved.

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To extract the tibia insert, the flat tip of the extractor is inserted into one of the slits left by the tibia insert upon being clipped to the tibial tray, rotating slightly (do not pry out).

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The patellar compression clamps, located on the upper tray of the patellofemoral set, is designed to make the bone cement brim up, and to maintain pressure during the hardening time necessary to properly attach the patella.

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The Genutech System has a tray impactor and a tibial component impactor. This latter instrument can impact the component once fully in position without damaging the functional areas of the tibial graft.

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
When positioning the femoral component, we can use the inserter/extractor for femoral components to guide the insertion and impact positioning of the component, then complete its fixation using the femoral impactor.

INTRODUCTOR / EXTRACTOR - INSERTER / EXTRACTOR

IMPACTOR FEMORAL / FEMORAL IMPACTOR

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To make minor adjustments in specific areas, we have a smaller impactor. This impactor, in conjunction with the first one, can be used to extract the test femoral components.

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There is a pin impactor which could be used as an extractor rod.

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*These components are manufactured on order. Contact Surgival.

Inserto tibial NPS
NPS tibial insert

	Talla Size	Longitud Length		Talla Size	Longitud Length
Ref. D8041100E	1	10 mm	Ref. D8044100E	4	10 mm
Ref. D8041120E	1	12 mm	Ref. D8044120E	4	12 mm
Ref. D8041140E	1	14 mm	Ref. D8044140E	4	14 mm
Ref. D8041160E	1	16 mm	Ref. D8044160E	4	16 mm
Ref. D8042100E	2	10 mm	Ref. D8045100E	5	10 mm
Ref. D8042120E	2	12 mm	Ref. D8045120E	5	12 mm
Ref. D8042140E	2	14 mm	Ref. D8045140E	5	14 mm
Ref. D8042160E	2	16 mm	Ref. D8045160E	5	16 mm
Ref. D8043100E	3	10 mm			
Ref. D8043120E	3	12 mm			
Ref. D8043140E	3	14 mm			
Ref. D8043160E	3	16 mm			

Implantes Implants

Inserto tibial PS
PS tibial insert

	Talla Size	Longitud Length		Talla Size	Longitud Length
Ref. D8051100E	1	10 mm	Ref. D8053160E	3	16 mm
Ref. D8051120E	1	12 mm	Ref. D8053180E	3	18 mm
Ref. D8051140E	1	14 mm	Ref. D8053200E	3	20 mm
Ref. D8051160E	1	16 mm	Ref. D8054100E	4	10 mm
Ref. D8051180E	1	18 mm	Ref. D8054120E	4	12 mm
Ref. D8051200E	1	20 mm	Ref. D8054140E	4	14 mm
Ref. D8052100E	2	10 mm	Ref. D8054160E	4	16 mm
Ref. D8052120E	2	12 mm	Ref. D8054180E	4	18 mm
Ref. D8052140E	2	14 mm	Ref. D8054200E	4	20 mm
Ref. D8052160E	2	16 mm	Ref. D8055100E	5	10 mm
Ref. D8052180E	2	18 mm	Ref. D8055120E	5	12 mm
Ref. D8052200E	2	20 mm	Ref. D8055140E	5	14 mm
Ref. D8053100E	3	10 mm	Ref. D8055160E	5	16 mm
Ref. D8053120E	3	12 mm	Ref. D8055180E	5	18 mm
Ref. D8053140E	3	14 mm	Ref. D8055200E	5	20 mm

Bandeja tibial
Tibial tray

	Talla Size
Ref. D8032100E	1
Ref. D8032200E	2
Ref. D8032300E	3
Ref. D8032400E	4
Ref. D8032500E	5

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Bandeja superior / Upper tray

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Longitud Length

Ref. D8120210	10 mm
Ref. D8120220	12 mm
Ref. D8120230	14 mm
Ref. D8120240	16 mm
Ref. D8120250	18 mm
Ref. D8120260	20 mm

4 Guía tibial simple
Simple tibial guide

Ref. D8220111

6 Soporte guía tibial (3 piezas)
Tibial guide support (3 pieces)


Ref. D8220110

8 Sistema fijación extramedular
Extramedullary fixation system

Ref. D8220150

9 Palpador tibial
Tibial probe

Ref. D8220270

2 Pin Ø3,4x55 mm con cabeza
Ø3,4x55 mm pin with head

Ref. D8210160

3 Pin Ø3,4x55 mm sin cabeza
Ø3,4x55 mm headless support pin

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
5 Cincel intramedular
Intramedullary chisel

Ref. D8220640

8 Destornillador 2,5 mm
2,5 mm screwdriver

Ref. D8220635


6 Impactor de quilla
Keel impactor

Ref. D8220460

9 Extractor de insertos
Inserts extractor

Ref. D8220620

Bandeja superior / Upper tray

1	Patela de prueba Trial patella	
		Diámetro Diameter
	Ref. D8130140	32 mm
	Ref. D8130150	34 mm
	Ref. D8130160	36 mm
	Ref. D8130170	38 mm
	Ref. D8130180	40 mm

2	Broca para patela Patellar drill bit
	Ref. D8230120

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2	Broca para tetones NPS NPS drill bit
	Ref. D8210230

5	Impactor femoral Femoral impactor
	Ref. D8210210

6	Pin Ø 3,4x40 mm con cabeza Ø 3,4x40 mm pin with head
	Ref. D8210162

4	Componentes femorales de prueba NPS NPS femoral trial components	
	Izquierdo Left	Talla Size Derecho Right
	Ref. D8111120	1 Ref. D8111110
	Ref. D8111220	2 Ref. D8111210
	Ref. D8111320	3 Ref. D8111310
	Ref. D8111420	4 Ref. D8111410
	Ref. D8111520	5 Ref. D8111510

D8601300 Set instrumental femoral 1 corte reducido / Reduced cutting femoral instrumental set 1


Bandeja superior / Upper tray

1 Destornillador 4,5 mm *
4,5 mm screwdriver
Ref. D8220630

6 Gancho palpador
Palpation hook
Ref. D8220290

8 Extractor de pins
Pin extractor
Ref. D8220280

2 Pin Ø 3,4x40 mm con cabeza
Ø 3,4x40 mm pin with head
Ref. D8210162

10 Medidor femoral primaria
Primary femoral sizing guide
Ref. D8271465

3 Pin Ø 3,4x55 mm con cabeza
Ø 3,4x55 mm pin with head
Ref. D8210160

4 Pinza extractora de pins
Pin extractor forceps
Ref. D8220285

11 Eje intramedular
Intramedullary shaft
Ref. D8210030

5 Punzón iniciador
Starter awl
Ref. D8210010

7 Broca de 9 mm
9 mm drill bit
Ref. 8210020

9 Barra auxiliar extractor/impactor de pins
Pin impactor/extractor auxiliar bar
Ref. D8271460

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
4 Base tibial
Tibial base plate
Ref. D8272510

7 Eje de alineación (400 mm)
Alignment shaft (400 mm)
Ref. D8210080

2 Pluma reducida primaria
Primary reduced stylus

Talla	Size
Ref. D8272215	1
Ref. D8272225	2
Ref. D8272235	3
Ref. D8272245	4
Ref. D8272255	5

9 Sistema de giro
Rotation system

Ref. D8272600	0°
Ref. D8272603	3°
Ref. D8272606	6°

3 Adaptador varo valgo
Varus-valgus adaptor


Ref. D8272335	3°
Ref. D8272355	5°
Ref. D8272375	7°

5 Acople adaptador rotación externo
Extractor rotation adaptor coupling
Ref. D8272340

8 Riel lateromedial
Latero-medial rail
Ref. D8272105


6 Mango roscado
Threaded handle
Ref. D8272520

D8601400 Set femoral instrumental 2 corte reducido / Reduced cutting femoral instrumental set 2



1	Máscara de corte biseles Bevel cutting mask		2	Introduccion/extractor componentes femorales Femoral component intramedullary/extractor		4	Impactor de componente femoral Femoral component impactor	
		Talla Size			Ref. D8271415			Ref. D8272530
	Ref. D8272051	1						
	Ref. D8272052	2						
	Ref. D8272053	3						
	Ref. D8272054	4						
3	Máscara de corte distal reducida primaria Reduced distal cutting mask primary		5	Máscara de recut distal reducida Reduced distal recutting mask				
		Ref. D8272050			Ref. D8272165			

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IR06EDT001 / 07-2015 © Rev. 1