KAREY®

Vástago KAREY®
KAREY® Stem

surgival
GRUPO COSÍAS

Características técnicas Technical details

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The cementless KAREY- HA primary stem is manufactured in titanium alloy (Ti6Al4V).

Straight stem of quadrangular section divided into two portions:

a) Proximal or metaphyseal, trapezoidal in shape on sagittal and frontal planes, obtaining three-dimensional primary stability and avoiding component subsidence.

b) Distal or diaphyseal, progressively subsized to reduce implant rigidity, improving contact with the cancellous bone, enhancing bone growth and secondary stability of the implant.

The progressive subsizing of the distal portion to avoid diaphyseal filling, reducing the "shielding stress".

The HA coating is applied using the latest technique of Air Plasma Spraying, with a thickness of 130 micron over the whole implant.

Both the HA coated cementless stem, and the KAREY-C cemented stem are implanted using the same instrumentation.

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The metaphyseal semicircular macrostructures increase the contact surface, avoiding component subsidence and having an antivaro effect.

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The longitudinal grooves increase the bone-implant contact surface, avoiding rotation and reducing distal rigidity with respect to the metaphyseal portion.

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Surgical approach • The objective of any surgical hip approach is correct visualization of the acetabulum and proximal femur. The choice of surgical hip approach will be based on the surgeon's preferences

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1 · Femoral neck osteotomy

With the femoral template on the preoperative X-ray, we will determine the proper femoral-neck level of resection. The osteotomy will be performed at a 45° dihedral angle (Fig. 1).

Fig. 2

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Fig. 4

2 · Preparation of the femur

Start opening the femoral canal using a hollow box chisel (Fig. 2), removing enough bone from the medial wall of the greater trochanter parallel to the femoral axis.

Insert the manual conical reamer (Fig. 3), as laterally as possible and parallel to the femoral axis to find the femoral canal.

With the smallest rasp mounted on the rasp-handle (Fig. 4), start the reaming, inserting it as far as the last line of teeth. This operation must be repeated increasing rasp size sequentially until the endomedullar cortical chatter is found. The anteversion is established with these rasping operations.

Fig. 5

Fig. 6

Fig. 7

Fig. 8

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3 · Trial reduction

With the rasp inserted in the femoral canal, carry out the reduction with the trial heads for rasps, in order to select the proper neck length (Fig. 5). Finally, once the neck length has been chosen, stability, tightness and leg length are tested.

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4 · Femoral component insertion

The definitive stem will correspond to the same size as the last rasp used. Mount the stem on the specific insertion handle (Fig. 6), applying impactation until it is inserted to the required depth, at this point the insertion handle is removed. With the final impactor the definitive stem can be impacted in a complementary way, if considered necessary (Fig. 7).

The trial femoral head for stem permit to do a new reduction to assure the length of the neck.

Finally we place the definitive head chosen on the stem and the articulation is reduced with the help of the impactor-reducer (Fig. 8).

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Note: The same instrumentation is used to implant both the KAREY HA stem and the KAREY C stem using the same technic shown. However, take into consideration that for the implantation of the KAREY-C one must consider the steps specific to all cemented stems: cavity cleaning, preparation and insertion of the bone cement, distal obturators, etc.


Instrumental Instruments

Set completo instrumental vástago KAREY
KAREY Stem instrumentation complete set

Ref. F0005005

Bandeja superior
Upper tray

1	Raspas KAREY KAREY Rasps
Ref. F0005209	Número 9 Number 9
Ref. F0005210	Número 10 Number 10
Ref. F0005211	Número 11 Number 11
Ref. F0005212	Número 12 Number 12
Ref. F0005213	Número 13 Number 13
Ref. F0005214	Número 14 Number 14
Ref. F0005215	Número 15 Number 15

2	Cabezas femorales de prueba para rasps Femoral trial heads for rasps
Ref. A1550040	Ø 28 C/C Short Neck
Ref. A1550041	Ø 28 C/M Medium Neck
Ref. A1550042	Ø 28 C/L Long Neck
Ref. A1550043	Ø 28 C/XL Extra Long Neck
Ref. A1550050	Ø 26 C/C Short Neck
Ref. A1550051	Ø 26 C/M Medium Neck
Ref. A1550052	Ø 26 C/L Long Neck
Ref. A1550053	Ø 26 C/XL Extra Long Neck

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Ref. A1536040	Ø 28 C/C Short Neck	Ref. A1536043	Ø 28 C/XL Extra Long Neck	Ref. A1536052	Ø 26 C/L Long Neck
Ref. A1536041	Ø 28 C/M Medium Neck	Ref. A1536050	Ø 26 C/C Short Neck	Ref. A1536053	Ø 26 C/XL Extra Long Neck
Ref. A1536042	Ø 28 C/L Long Neck	Ref. A1536051	Ø 26 C/M Medium Neck		

Bandeja inferior
Lower tray

1	Impactor reductor Impactor-reducer
Ref. A2700375	
2	Fresa iniciadora endomedular Endomedullary starting reamer
Ref. A2700600	
3	Cinzel hueco nº 2 Hollow box chisel nº 2
Ref. A2300805	
4	Extractor para cabeza femoral Femoral head extractor
Ref. A1150010	

5	Punzón impactador para vástago Impactor awl for stem
Ref. F0005500	
6	Mango impactor vástago definitivo Impactor handle definitive stem
Ref. F0005300	
7	Varilla para mango de raspa Tommy bar for rasp handle
Ref. A1530101	
8	Mango automático para rasps Automatic handles for rasps
Ref. A1555005	
2 unidades 2 units	

Implantes Implants


Karey-HA

Num. 9	Ref. F0005009E
Num. 10	Ref. F0005010E
Num. 11	Ref. F0005011E
Num. 12	Ref. F0005012E
Num. 13	Ref. F0005013E
Num. 14	Ref. F0005014E
Num. 15	Ref. F0005015E

Karey-C

Num. 9	Ref. F0005109E
Num. 10	Ref. F0005110E
Num. 11	Ref. F0005111E
Num. 12	Ref. F0005112E
Num. 13	Ref. F0005113E
Num. 14	Ref. F0005114E
Num. 15	Ref. F0005115E

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Cabezas femorales de acero inoxidable*
Stainless-steel femoral heads*

Ref. A1507040E	Ø 28 C/C Short Neck
Ref. A1507041E	Ø 28 C/M Medium Neck
Ref. A1507042E	Ø 28 C/L Long Neck
Ref. A1507050E	Ø 26 C/C Short Neck
Ref. A1507051E	Ø 26 C/M Medium Neck
Ref. A1507052E	Ø 26 C/L Long Neck
Ref. A1509040E	Ø 28 C/C Short Neck
Ref. A1509041E	Ø 28 C/M Medium Neck
Ref. A1509042E	Ø 28 C/L Long Neck
Ref. A1509043E	Ø 28 C/XL Extra Long Neck
Ref. A1509050E	Ø 26 C/C Short Neck
Ref. A1509051E	Ø 26 C/M Medium Neck
Ref. A1509052E	Ø 26 C/L Long Neck
Ref. A1509053E	Ø 26 C/XL Extra Long Neck

*Disponible cabeza femoral en Cromo Cobalto
*Avaliable femoral head Cobalt-Chrome

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ICAKACT01 / 03-2015 © Rev. 1