*Técnica Quirúrgica
Surgical Technique*

SHERPA®

*Sistema de Columna SHERPA®
SHERPA® Spine System*

 **surgival**
GRUPO COSÍAS

Índice

Table of contents

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It is advisable to make a central posterior craniocaudal incision, which can be made longer or shorter depending on the levels to be fused.

It is important to skeletonise the vertebral segments that require intervention very carefully, minimising bleeding in the area in order to identify the different anatomical levels easily.

The segments of the posterior arch are dissected from the medial line to the tip of the transverse process.

The technique used to prepare the pedicle differs depending on whether you are working with sacral, lumbar or thoracic vertebra.

When locating the point of entry for the screw into the pedicle, the image amplifier should be used.

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Initially, the superior and transverse articular processes on the same level must be skeletonised. The entry point on the pedicle usually coincides with the intersection between the medial line of the transverse process and the base of the articular process.

To prevent the instrument from slipping on the surface of the articular process, it is advisable to remove the bone cortex from this area using rongeur forceps, a chisel or motorised reamer.

Particular care must be taken with the immediately superior vertebra so as not to damage its inferior articular process. Then, working in a direction convergent with the medial line and at an angle of around 30°, the path required for the screw must be hollowed out using the specific instruments provided for this purpose: awl, finder, feeler, tap and curette.


***NOTE:** When using the pedicle finder, the first 20 mm can be inserted with the point towards the lateral cortex of the pedicle; it can then be removed and reinserted at an angle of 180° towards the internal cortex to proceed with the preparation.

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A small osteotomy is made, oblique to the distal end of the inferior articular process of L5, using a chisel, to visualise the articular process of the sacrum. Saucerisation or spongialisation may be performed using rongeur forceps or a chisel. At the base of the articular process there is a small bone crest or accessory process, which is considered the best point of entry for the sacral pedicle. The awl should be used to facilitate its location.

Then, working in a direction convergent to the medial line, the path to be followed by the S1 pedicle screw is hollowed out using the curette, pedicle finder and the tap feeler. The best direction for the sacral screw is towards the most anterior and superior part of the S1 vertebral endplate.

Using the pedicle feeler, you can feel the end and the walls of the path created.

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Note: The same instruments are used as for the lumbar pedicle

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The posterior arch must be exposed sufficiently to the distal end of the transverse process. The point of entry of the pedicle in the thoracic vertebra coincides with the join between the lines traced in the upper third of the transverse process and the outermost third of the thoracic plate on the medial line of the articular process.

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The technique used to expose the area is similar to that described above. To access the entry point, the bone cortex posterior to the transverse process should be removed (with rongeur forceps, for example), between the tip and the middle third.

The curette is inserted through the hole made, in a convergent direction at an angle of around 20°, in the same direction as the transverse process, until you penetrate to the vertebral body in the area anterior to the pedicle, traversing the costotransverse joint.

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Note: The same instruments are used as for the lumbar pedicle.

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When inserting the transpedicular screw, having previously chosen the required length and diameter, the following factors must be taken into consideration:

- The diameter of the pedicle increases in the craniocaudal direction (L1-L5).
- The angle of convergence of the lumbar pedicle increases from L1 through to L5 (between 5 and 20 degrees).
- The chosen length must not perforate the anterior bone cortex of the vertebral body.
- For most thoracic segments, the pedicle has a considerable height but a small transverse diameter (for many levels, less than 5 mm).
- The extrapedicular trajectory is recommended for high thoracic segments, from T9 - T8 to T2. Bear in mind that the outer ring of the screw will be more lateralised than in the case of the pedicular trajectory.

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Having prepared and checked the trajectory of the pedicle using the feeler, the screw is fitted into its specific key. The proximal end of the key is inserted into the slots on the outer ring and the external component is threaded into the internal surface of the ring. The screw is placed in its position and, finally, the insertion key is removed by unscrewing it from the outer ring.

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

Having located the path of the pedicle, the screw is fitted into the specific screwdriver. The internal part of the key fits into the hexagon on the head, and the external part threads into the internal surface of the outer ring. The screw is placed in its position, and then the insertion key is removed by unscrewing it from the outer ring.

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

Using the rod holders, offer the rod up to the screws to select the right length. At the same time, you will get an image of the outline required to start moulding the rod. The rod should be moulded using the bending pliers, in position 6, equivalent to the Ø 6 mm rod, and the rod holders.

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Next, make sure the rod is sitting right at the back of all the outer rings, using the rod pusher.

If the screws are positioned so close to one another that the outer rings are touching, you may use another pusher inside the outer rings.


If you cannot achieve a snug fit between the rod and the back of the outer ring with the pusher, you can use the outer ring lever.

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Having inserted the rod, lock the system in place using the various tightening elements, depending on whether you are locking a monoaxial, polyaxial or spondylolisthesis screw.

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To fix the rod in the monoaxial screw, use the bevelled washer and the monoaxial system cap.

The bevelled washer slides onto the washer inserter, which has been previously threaded onto the screw.

Then, slide the washer guide onto the same inserter. This instrument is used to guide the washer, by means of rotation, so that it rests on the cranial and caudal ends of the rod.

Finally, insert the cap, which will fix the system in place.

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
NOTE: Use the same precautions as polyaxial screw locking (see next page)

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

Fix the rod to the polyaxial screw using the polyaxial cap.

The polyaxial cap is threaded onto the polyaxial cap to transport it to the outer ring where it is threaded inside until it lies flush to the ring.

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Use the spondylolisthesis cap and the external washer to fix the bar to the spondylolisthesis screw.

Without removing the polyaxial cap rod, place the external washer onto the spondylolisthesis washer key and slide it down the key.

The external washer will be in the right position to be threaded onto the outside of the outer ring, alternately inserting the cap and the external washer until contact is made with the rod.

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When tightening the locking systems on the rod, the outer ring must be placed parallel to the spinous process on the sagittal plane so that, in traction, there is only anteroposterior displacement.

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Once retraction is complete, the outer rim must be split using the outer rim splitter. Make sure that both the washer and the screw are fully screwed in place.

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Compression or distraction may be performed using the specific forceps provided for these tasks before carrying out the final tightening stage.

Pinza compresora
Compression forceps

Pinza distractora
Distractor forceps

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To prevent the screws from becoming unfastened when they are tightened, the anti-turn key may be used.

Llave antigiro
Anti-turn key

Destornillador hexagonal
Hexagonal screwdriver

Implantes Implants

Set implante Sherpa®
Sherpa® implant set

Ref. C6836000


Barra Sherpa® Ø 6 mm Rod Sherpa® Ø 6 mm	Longitud Length
Ref. C6800740	40 mm
Ref. C6800750	50 mm
Ref. C6800760	60 mm
Ref. C6800770	70 mm
Ref. C6800780	80 mm
Ref. C6800790	90 mm
Ref. C6800800	100 mm
Ref. C6800810	110 mm
Ref. C6800820	120 mm
Ref. C6800830	130 mm
Ref. C6800840	140 mm
Ref. C6800850	150 mm
Ref. C6800860	160 mm
Ref. C6800870	170 mm
Ref. C6800880	180 mm

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ICOSHET001/ 07-2015 © Rev. 1