SHINE-C

Vástago Cementado Shine-C
Shine-C Cemented Stem

surgival
GRUPO COSÍAS

Características Técnicas *Technical Details*

- Cono 12/14
- Offset constante 41 mm.
- Eje cervicodiafisario de 130°
- 12/14 morse neck
- Unique offset 41 mm.
- 130° neck-shaft angle

Cotilos cementados *Cemented cups*

Con ceja de 15°
Antiluxation offset 15°

Standard plano
Standard flat cup

· Cabeza inox. 28 mm. Ø
· 28 mm. Ø Stainless-steel head

· Acero nitrogenado
· Nitrogen-steel

· C-Plug® Restrictor de cemento*
· C-Plug® Cement restrictor*

*Consultar referencias con el distribuidor o fabricante
*Consult your distributor or manufacturer for references

Artroplastia primaria cementada de cadera

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1 · Patient positioning

The surgical procedure can be done with the patient in the lateral position, or in supine decubitus. The lateral position is better suited for the posterolateral approach, for obese patients, for revision procedures and for difficult primary procedures in which trochanteric osteotomy will be needed.

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2 · Femoral neck osteotomy (Fig. 1)

With the femoral template on the preoperative X-ray, we will know the proper femoral-neck level of resection by measuring this distance from the lesser trochanter.

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3 · Acetabular exposure and reaming (Fig. 2)

A full, and wide view of the whole acetabulum, is essential to prepare it in the right way (removing the soft tissues, reaming the cartilage, drilling the holes into the cancellous bone), and to place correctly the acetabular cup. By placing Hohmann's retractors over the anterior pelvic brim, under the posterior acetabular wall and below the acetabular notch, we get a full view of the acetabulum.

Acetabular reamers in increasing size, are used to remove all the cartilage tissue and to obtain a raw bleeding bone.

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4 · Acetabular drilling

A stop-drill is then used to make several holes into the acetabular bone. These holes will later be filled with acrylic cement.

Cemented primary hip arthroplasty

The SHINE-C is a collarless and highly polished surface femoral stem. It is nowadays very well accepted that the collarless and polished stems minimizes the risk of stem loosening and of wear debris production. This type of stem is one of the first options in cases of primary hip arthroplasties.

Accurate preoperative planning is mandatory in order to obtain a successful result. By using the femoral templates provided by the manufacturer, the surgeon will know:

- The femoral stem size.
- The level of the neck osteotomy.
- The appropriate neck-length of the femoral head.

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5 · Cementation of the acetabular cup (Fig. 3)

After thoroughly cleaning and rinsing the acetabular cavity, a portion of bone cement in doughy state is forcefully introduced into the acetabular cavity. The polyethylene cup is mounted into the acetabular introductory guide. The polyethylene cup is placed into the acetabular cavity with the appropriate inclination and anteversion angle. The acetabular cup guide is pushed forward to seat the cup against the acetabular bone wall and to remove the extra cement. The construct is kept in place until the cement is completely set.

The final polyethylene cup must be 2-3 mm. narrower than the acetabular cavity to get a proper cement wall thickness around it.

Fig. 4

Fig. 5

Fig. 6

Fig. 7

Fig. 8

Fig. 9

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6 · Preparation of the femur

A. A Box Osteotome is hammered into the femoral neck osteotomy, parallel to the femoral canal (Fig. 4). Enough amount of the medial wall of the greater trochanter to reach its tip, is removed.

B. After removing this piece of bone, a straight and long hand-driven awl or reamer, is introduced as far lateral as possible, into the femoral canal to identify it (Fig 5).

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7 · Rasp insertion

A. The SHINE-C rasps are mounted on the rasp-handle (Fig. 6). It is advisable to start the rasping with one or two sizes smaller than the size selected during the templating. The rasp is kept parallel to the anterior and posterior neck cortices, to reproduce the previous neck anteversion. The size of the rasps are sequentially increased until cortical chatter is found. The upper end of the rasps show three visible lines to know the depth of insertion. Once the proper depth is obtained, we are ready to make the first trial reduction.

B. Trial heads with different neck lengths are fitted over the rasp neck (Fig. 7). A trial reduction is made. Joint stability, muscle tightness and leg length are tested.

If we are satisfied with the head-neck-length selected, the rasp is removed and the femoral canal is prepared for cementing. It is recommended to plug distally the femoral canal with a cement restrictor to pressurize the cement. This plug is located 2 cm. distally to the tip of the femoral stem. The cement can be introduced into the femoral canal by hand or using a cement gun.

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8 · Femoral component insertion

The femoral component of the selected size is mounted over its handle (Fig. 8). The femoral stem is manually introduced into the cemented femoral canal until the desired depth.

When the bone cement is set, we can make a second head-neck size trial reduction by fitting the different trial head-necks over the neck of the cemented stem. (Fig. 9).

If we are satisfied with the head-neck length chosen, the final head is sharply tapped over the stem neck. The hip joint is reduced.

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
Ref. A1536040	Ø 28 C/C	Short Neck
Ref. A1536041	Ø 28 C/M	Medium Neck
Ref. A1536042	Ø 28 C/L	Long Neck
Ref. A1536043	Ø 28 C/XXL	Extra Long Neck
Ref. A1536050	Ø 26 C/C	Short Neck
Ref. A1536051	Ø 26 C/M	Medium Neck
Ref. A1536052	Ø 26 C/L	Long Neck
Ref. A1536053	Ø 26 C/XL	Extra Long Neck

Bandeja superior Upper tray

2 Cabezales femorales de prueba para rasps Femoral trial heads for rasps

Ref. A1550040	Ø 28 C/C	Short Neck
Ref. A1550041	Ø 28 C/M	Medium Neck
Ref. A1550042	Ø 28 C/L	Long Neck
Ref. A1550043	Ø 28 C/XXL	Extra Long Neck
Ref. A1550050	Ø 26 C/C	Short Neck
Ref. A1550051	Ø 26 C/M	Medium Neck
Ref. A1550052	Ø 26 C/L	Long Neck
Ref. A1550053	Ø 26 C/XL	Extra Long Neck

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
Cabezales femorales acero inoxidable Stainless-steel femoral heads

Ref. A1509040E	Ø 28 C/C	Short Neck
Ref. A1509041E	Ø 28 C/M	Medium Neck
Ref. A1509042E	Ø 28 C/L	Long Neck
Ref. A1509043E	Ø 28 C/XL	Extra Long Neck

Cotillos cementados Cemented cups

	15°	mm. Ø	Standard
Ref. A1512244E	28 x 44		Ref. A1511244E
Ref. A1512246E	28 x 46		Ref. A1511246E
Ref. A1512248E	28 x 48		Ref. A1511248E
Ref. A1512250E	28 x 50		Ref. A1511250E
Ref. A1512252E	28 x 52		Ref. A1511252E
Ref. A1512254E	28 x 54		Ref. A1511254E
Ref. A1512256E	28 x 56		Ref. A1511256E
Ref. A1512258E	28 x 58		Ref. A1511258E

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