SHY®

Cotilo S.H.Y.®
S.H.Y.® Cup

surgival
GRUPO COSÍAS

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The Modular Acetabular System or S.H.Y.® Cup from Surgival is made up of the following subcomponents:

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The hemispheric cup is made from titanium (Ti6Al4V). Primary attachment is provided by the ridged section around the middle, for more efficient attachment, and Ø 6.5 mm low profile titanium screws.

For secondary attachment, this cup is coated in titanium plasma spray to increase the roughness and contact surface of the dome. Finally, it is covered with a complete coating of Hydroxyapatite (HA), to foster the osteoinduction and osseointegration of the bone on the treated surface.

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
1. The top hole provides a viewing point to ascertain contact between the implant and the bone.
2. The three holes enable additional attachments to be made to the primary attachment using Ø 6.5 mm low profile titanium screws, to be used when the surgeon determines that an acceptable level of 'pressfit' stability has not been achieved.
3. The peripheral slots around the base perform various functions:
 - a) Provide primary 'pressfit' stability by increasing the diameter of the base.
 - b) Increase the bone-implant surface area treated with Hydroxyapatite to foster osteoinduction and osseointegration, providing optimum biomechanical anchoring.
 - c) Decrease micro-movements of the implant on the bone surface.
4. The extraction slots around the base facilitate the removal of the Insert.
5. Internal diameter groove to anchor the implant.

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6. The insert should not come above the line of the cup, with the exception of the 15° anti-luxation lip, in order to maximize the possible range of motion of the prosthesis. The width of the Insert will always be over 5 mm, even on smaller sizes.

7. System designed with a conical trunk section offering maximum congruence, allowing for a greater spread of tension between insert/cup interpaces. This will increase the durability of the implant by reducing the stresses to which polyethylene is usually subjected.

8. The Insert is 'clipped' to the Cup and a 'harpoon' type flange is used to prevent it from becoming loose. The system clips onto the diameter groove of the metal cup.

9. It has 12 peripheral notches, spaced 30° apart, which control rotation, resist inclination and luxation forces, and offer greater versatility during insert placement.

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
10. Self-threading tip, with 3 cutting edges.

11. Maximum angling of 15° so the implant can be attached to areas of the bone offering a greater structural quality.

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1 · Patient's positioning and surgical approach

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The patient can be placed either in the lateral or in the supine decubitus, according to the surgeon's preferences. The lateral decubitus is most frequently used in all the cases where a posterolateral approach should be done. For the above approaches, anterolateral and lateral approaches can be chosen for the supine or lateral decubitus positions.

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With the help of Hohmann's retractors or Steinmann nails placed in the acetabular lips, a wide exposure of the acetabulum should be obtained.

Cleaning all the soft tissues from inside the acetabular cavity, including the synovial membrane, the ligamentum teres and the inferior transverse ligament, allows a complete and wide exposure of the acetabulum.

By using increasing size reamers, the articular cartilage and a thin layer of subchondral bone is removed, until a regular bleeding surface is obtained. The least possible amount of bone will be milled to maintain a good bone stock and not weaken the acetabular wall, but it should have complete coverage by the reamer. The position of the milling is at 45° of lateral inclination and 15° of anteversion.

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In most of the cases the chosen acetabular cup matches the size of the last acetabular reamer. The acetabular cup is oversized by two mm on the most peripheral part, to allow a stable fit after forceful fixation.

The trial acetabular components give very useful information about the size (diameter) and the quality of the inner and superolateral fitting of the acetabular component which is best adapted to each patient. Through its windows we can see the degree of adaptation to the base of the acetabular cavity.

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The acetabular cup is assembled on the acetabular guide according to the surgeon's criteria and is introduced with the proper 45° of inclination and 15° of anteversion.

Sometimes, after impaction of the cup, the surgeon wants to improve the degree of stability of the component. Screws can be threaded with a range of 15° angulations through the component and into the acetabular bone.

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Position of the guide rods:

The rod that controls 45°, if the patient is in the lateral decubitus position, must be perpendicular to the floor. If the patient is in the supine decubitus position, the rod that controls the 45° must remain perpendicular to the longitudinal axis of the patient. In this position, 15° of anteversion will be given to the impactor.

(When in doubt, please request the Cup insertion's guide).

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The polyethylene cup can be placed with its anti dislocation lip in a wide range of different 30° angulations.

The different positions are obtained by matching the peripheral geometry of the polyethylene and metallic cups.

The removal of the polyethylene component is done by using an extractor and levering through a superior notch in the metallic cup.

The anti-dislocation lip is used to avoid dislocations. In the cases in which a posterior approach has been used, if dislocation occurs it will most probably be posterior. For this reason, in this situation, it is advisable to place the lip in the posterosuperior position.

In the cases in which an anterior or anterolateral approach has been used, the dislocation will preferably occur on the anterior path.

Therefore, for this situation, it would be advisable to place the lip at 45° and in a neutral position of anteversion-retroversion.

The anti-dislocation lip in anterior position is not usually advisable because it can produce interferences with the flexion movements of the prosthesis, by colliding with the femoral neck.

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*For the removal of the polyethylene insert, standard extractors will be used. This may be difficult due to the adjustment between the insert and the acetabular component.

It is possible that the polyethylene insert may be damaged during its removal. In this case it must be replaced by a new one.

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For the correct assembly of the impactor, the tabs (indicated by two arrows in the lower part, fig. 1a) must fit perfectly within the slots of the acetabular component (fig. 1b) until verifying that both component and impactor turn at the same time (fig. 2).

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Next, press down the on the wheel located in the upper part (fig. 3) and turn it so that the threads of the impactor are completely inserted in the acetabular component. Once this step is done and the fastening is ensured, the guide rods (fig. 4) are installed and the impactation is carried out.

Nota: si el montaje no se realiza exactamente como se ha indicado, se corre el riesgo de que las pestañas o el tornillo superior puedan llegar a dañarse, inhabilitando por completo el funcionamiento del impactor.
Note: if the assembly is not done exactly as indicated, there is a risk that the tabs or the upper screw may become damaged, completely disabling the functioning of the impactor.

Implantes Implants

Tornillos Screws

Ref. A2400520	20 mm longitud	length
Ref. A2400525	25 mm longitud	length
Ref. A2400530	30 mm longitud	length
Ref. A2400535	35 mm longitud	length
Ref. A2400540	40 mm longitud	length

Insertos cotilo Cup inserts

Ref. A2400746E	Ø 46 mm / Ø 26 mm
Ref. A2400748E	Ø 48 mm / Ø 26 mm
Ref. A2400750E	Ø 50 mm / Ø 28 mm
Ref. A2400752E	Ø 52 mm / Ø 28 mm
Ref. A2400754E	Ø 54 mm / Ø 28 mm
Ref. A2400756E	Ø 56 mm / Ø 28 mm
Ref. A2400758E	Ø 58 mm / Ø 28 mm
Ref. A2400760E	Ø 60 mm / Ø 28 mm

Cotilos Cups

Ref. A2400646E	Ø 46 mm HA
Ref. A2400648E	Ø 48 mm HA
Ref. A2400650E	Ø 50 mm HA
Ref. A2400652E	Ø 52 mm HA
Ref. A2400654E	Ø 54 mm HA
Ref. A2400656E	Ø 56 mm HA
Ref. A2400658E	Ø 58 mm HA
Ref. A2400660E	Ø 60 mm HA

Set completo instrumental S.H.Y.* S.H.Y. instrumentation complete set*

Ref. A2400001

Bandeja inferior Lower tray

1	Cotilos de prueba Trial acetabular cups	
Ref. F0005846	Ø 46 mm	
Ref. F0005848	Ø 48 mm	
Ref. F0005850	Ø 50 mm	
Ref. F0005852	Ø 52 mm	
Ref. F0005854	Ø 54 mm	
Ref. F0005856	Ø 56 mm	
Ref. F0005858	Ø 58 mm	
Ref. F0005860	Ø 60 mm	
Ref. F0005862	Ø 62 mm	


2	Impactor definitivo inserto-cotilo Definitive cup-insert impactor	Ref. A2400834	Ø 26 mm
		Ref. A2400835	Ø 28 mm
3	Mango para cotilo Cup's handle	Ref. F0005835	
4	Extractor para inserto Inserts extractor	Ref. F0005900	
	2 unidades 2 units		

5	Posicionador inserto Positioner for inserts	
Ref. A2400831	46/48 Ø 26 mm	
Ref. A2400832	50/52 Ø 28 mm	
Ref. A2400833	54/62 Ø 28 mm	

6	Impactor orientador cotilo Cup guide impactor	Ref. A2400830
7	Impactor final para cotilo Cup final impactor	Ref. A2400870

Set instrumental Cotilo S.H.Y. S.H.Y. cup instrumentation set


Ref. A2400005

Bandejas Cotilo S.H.Y. S.H.Y. Cup trays

Ref. A2400940

Bandeja superior Upper tray

1	Destornilladores hexagonales Hexagonal screwdrivers	
Ref. B3610000	Recto Straight	
Ref. B3611350	Cardan Cardan	
2	Pinza porta-tornillos Screw holding forceps	Ref. F0005950
3	Guía de broca Drill guide	Ref. F0005935

4	Árbol flexible para brocas Drills flexible shafts	Ref. F0005930
	2 unidades 2 units	
5	Brocas Drill bits	
Ref. F0005931	Ø 3,2 x 35 mm	
Ref. F0005932	Ø 3,2 x 56 mm	
6	Medidor de profundidad Depth gauge	Ref. F0005960

Set fresas acetabulares Acetabular reamers set

Ref. A1701000

1	Fresa acetabular Acetabular reamer	
Ref. A1701044	Ø 44 mm	
Ref. A1701046	Ø 46 mm	
Ref. A1701048	Ø 48 mm	
Ref. A1701050	Ø 50 mm	
Ref. A1701052	Ø 52 mm	
Ref. A1701054	Ø 54 mm	

Ref. A1701056	Ø 56 mm	
Ref. A1701058	Ø 58 mm	
Ref. A1701060	Ø 60 mm	
Ref. A1701062	Ø 62 mm	
2	Acople universal Universal Join	Ref. A1700380
3	Mango porta fresa acetabular universal Handle for universal acetabular reamer	
Ref. A1701075		
	2 unidades 2 units	


Caja set fresas acetabulares cotilo S.H.Y. S.H.Y. Cup box for acetabular reamer set

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ICASHVCT01 / 07-2015 © Rev. 1