
SURGIFEET®

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DIAPHYSEAL OSTEOTOMY OF THE 1ST METATARSAL. **“SCARF” TECHNIQUE**

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Internal access and exposure of the metatarsophalangeal joint.

Carry out the lateral release of the metatarsal-phalangeal-sesamoid complex. This release can be carried out via the first dorsal commissure or through the internal access. Cut the suspensory ligament of the sesamoid bones. This will allow the sesamoid bones to move towards the sole, being placed on the head of the metatarsal.

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a) Internal incision at the junction between the dorsal and plantar skin, until the bone plane and joint release.

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b) Preserve the lateral dorsal capsule to maintain the bloody supply via the capital dorsal branch of the metaphysis and the plantar capsule by not performing the incision at the entrance of the distal plantar vessels.

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
c) Exostectomy outside the sagittal groove of the metatarsal head.

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
2 Kirschner wires Ø 0.9 mm are inserted.

- The first wire is inserted in the distal dorsal position of the 1st metatarsal angled slightly towards the plantar position.
- The second wire is inserted in the proximal position, in the metaphyso-diaphyseal area parallel to the 1st wire, and they may be slightly convergent. These 2 wires show the level at which the osteotomy will be performed.

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
The diaphyseal line must end near to the plantar side of the proximal end of the bone (always below the medial crest, which serves as a point of reference) and near to the dorsal side of its distal end.

The correct position of these sections of the line of the osteotomy will increase its strength during dorsiflexion.

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Using these wires as a reference point, cut lengthwise along the diaphysis, angled from dorsal-distal to plantar-proximal.

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
Cut the wires to prevent problems of accessibility with the motor. When cutting the diaphysis, you can use the depth gauge to help guide the cutting.

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The distal cut is made at an angle of 60° to the plane of the diaphyseal osteotomy, following the direction of the Kirschner wire which has been placed there as a point of reference. The proximal cut is parallel to the distal cut and follows the direction of the proximal Kirschner wire.

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
If it is necessary to cut the 1st metatarsal there are 2 possibilities:

- Make the transverse cuts on an angle.
- Resect the ends of the proximal and distal fragments.

Shortening the 1st metatarsal decompresses the metatarsophalangeal joint.

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
Then perform a lateral displacement and correct the proximal articular set angle (PASA). For each mm of lateral displacement there is a 1° correction of the intermetatarsal angle.

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
Now insert a Kirschner wire with a diameter of Ø 0.9 mm.

The point of entry for the 1st screw is located 1 cm away from the osteotomy and on the medio-lateral plane; the point of insertion is at the junction between the internal 2/3 and the external 1/3 of the diaphysis.


The 2nd screw is inserted in the midpoint between the entry point of the 1st screw and the proximal line of the osteotomy. On the medio-lateral plane, the point of entry must be located in the middle of the diaphysis. Thus, strength will be increased for dorsiflexion originating in the proximal area of the metatarsal.

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
Insert the Kirschner wire, aiming it at the centre of the head of the 1st metatarsal until it can be seen to exit through the distal plantar metatarsal area, remaining flush in the head of the screw. Then take the measurement with the depth gauge.

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The screw to be used must be 7 mm shorter than measurement obtained to avoid damaging the cartilage and the metatarsophalangeal joint.

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
Then use the drill for the core of the SURGIFEET® screw. This drill is cannulated and has a depth gauge system. Therefore, drill the length of the selected screw.

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If the drill is the same length as the result obtained from the gauge there is a risk of losing the Kirschner wire as a point of reference for the subsequent steps and damaging the articular cartilage.

Tras este primer fresado, se realiza un escariado para el núcleo de la cabeza del tornillo.
Following this initial drilling, perform reaming for the core of the head of the screw.

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Then, insert the distal screw. When selecting the screw take into account the foot which is undergoing surgery.

The SURGIFEET® system has screws which are right-hand threaded (dextrorotatory) and left-hand threaded (levorotatory):

- The right-hand threaded (dextrorotatory) screws are used for left feet.
- The left-hand threaded (levorotatory) screws are used for right feet.

This system prevents undesired rotation of the “SCARF” metatarsal osteotomy when the screws are tightened.

The distal screw is inserted with a screwdriver with a tapered hexagonal fastening system. The screw is taken out of its instrument cabinet with this screwdriver and it remains preassembled until it is inserted.

Insert until it is flush or slightly inside the bone.

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Insert the Kirschner wire obliquely until it can be seen to exit through the external plantar metatarsal area. Then take the measurement with the depth gauge.

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The screw to be selected must be 2 mm shorter than measurement obtained.

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
Then select the drill for the core of the SURGIFEET® screw. This drill is cannulated and has a depth gauge system.

Before drilling, insert the Kirschner wire 4 mm above the plantar area. Then clamp it with mosquito forceps to prevent it from moving during subsequent drilling and insert the screw.

Drill the cortices of the proximal and distal fragments with the drill.


Tras este primer fresado, se realiza un escariado para el núcleo de la cabeza de tornillo. Finalmente, se introduce el tornillo seleccionado.
Following this initial drilling, perform reaming for the core of the head of the screw. Finally, insert the selected screw.

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Then resect the surplus bone fragment of the diaphysis. Additionally, check the stability and correction obtained, both on the dorsoplantar plane and the lateral plane.
Finally, release the median nerve and perform a capsulorrhaphy.

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DISTAL OSTEOTOMY OF THE 1ST METATARSAL. “CHEVRON” TECHNIQUE

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
This is performed in the same way as the “SCARF” technique, bearing in mind that the incision in the skin reaches from the distal part of the metatarsal to the middle of the proximal phalanx, and it is strictly median so as to not damage the digital nerve.

Make a longitudinal incision in the capsule until you reach the bone plane and joint plane. Expose the metatarsal head. During this step you must take care not to excessively raise the periosteum of the metatarsal head and you must preserve plantar vascularisation.

The median eminence or exostosis is resected with an oscillating saw or fine osteotomy parallel to the diaphyseal axis and just at the level of the groove which delimits it.

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Locate the point of entry of the Kirschner wire 2 mm from the dorsal area and 4 mm proximal to the articular surface. The Kirschner wire must be perpendicular to the 2nd metatarsal and angled slightly towards the sole of the foot.

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If it is necessary to cut the 1st metatarsal there are 2 possibilities:

- a) Make the transverse cuts on an angle.
- b) Resection of the proximal end.


Shortening the 1st metatarsal decompresses the metatarsophalangeal joint.


If the PASA needs to be corrected you can perform a wedge resection with an internal base.

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An L-shaped osteotomy is performed with a thin, narrow and straight saw. The proximal line of the "L" must follow the plane defined by the Kirschner wire, joining the point of insertion of the wire with the proximomedial area of the diaphysis. The distal line is made 60° from the proximal line, following the direction of the reference wire.

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

When the osteotomy has been completed, perform the lateral displacement of the head. For every millimetre of displacement you can expect a 1° decrease in the IMA (Intermetatarsal Angle). Therefore, this osteotomy is indicated for patients with an IMA below 14°. It is desirable for the displacement to range from 4 to 6 mm at most.

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
Now insert a Kirschner wire with a diameter of Ø 0.9 mm. The point of entry for the screw is located 1 cm away from the osteotomy and on the medio-lateral plane; the point of insertion is at the junction between the internal 2/3 and the external 1/3 of the diaphysis.

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Insert the Kirschner wire, aiming it at the centre of the head of the 1st metatarsal until it can be seen to exit through the distal plantar metatarsal area, remaining flush in the head of the screw.

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Then take the measurement with the depth gauge.

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If the drill is the same length as the result obtained from the gauge there is a risk of losing the Kirschner wire as a point of reference for the subsequent steps and damaging the cartilage of the head.

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
Following this initial drilling, perform reaming for the core of the head of the screw.

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The diameter of the head is Ø 3.35 mm and the diameter of the screw is Ø 2.55 mm.

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Then, insert the distal screw. When selecting the screw take into account the foot which is undergoing surgery. The SURGIFEET® system has screws which are right-hand threaded (dextrorotatory) and left-hand threaded (levorotatory):

- The right-hand threaded (dextrorotatory) screws are used for left feet.
- The left-hand threaded (levorotatory) screws are used for right feet.

This system prevents undesired rotation of the “CHEVRON” metatarsal osteotomy when the screws are tightened.


The distal screw is inserted with a screwdriver with a tapered hexagonal fastening system. The screw is taken out of its instrument cabinet with this screwdriver and it remains preassembled until it is inserted, until flush or slightly inside the bone.

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Then resect the surplus bone fragment of the diaphysis.

Additionally, check the stability and correction obtained, both on the dorsoplantar plane and the lateral plane.

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It is recommended to make a median incision of around 6 cm between the dorsal and plantar skin, starting at the level of the interphalangeal joint. Slightly raise the periosteum of the first phalanx. At this point you must have located the tendon of the flexor hallucis longus to avoid damaging it.

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
Perform an osteotomy which is perpendicular to the axis of the proximal phalanx at the metaphyseal-diaphyseal junction.

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
Then perform a second osteotomy which is distal to the first with the angulation you want to correct. The external cortex must remain intact so that it can perform hinge-like movements.

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
When the osteotomy has been performed, select which type of staple to use, straight or oblique. Mount it on the staple holder-impactor and insert it.

Finalmente, se dispone de un impactador final de grapas con el que impactar la misma hasta su posición final.
Finally, there is a final staple impactor with which it is impacted until it reaches its final position.

Implantes Implants

Tornillo canulado para pie izquierdo
(Rosca a derecha)
Cannulated screw for left foot
(Right-hand threaded)

	Longitud Length	Unidades Units
Ref. B3324010	10 mm	3
Ref. B3324012	12 mm	3
Ref. B3324014	14 mm	3
Ref. B3324016	16 mm	7
Ref. B3324018	18 mm	7
Ref. B3324020	20 mm	7
Ref. B3324022	22 mm	3
Ref. B3324024	24 mm	3
Ref. B3324026	26 mm	3
Ref. B3324028	28 mm	3
Ref. B3324030	30 mm	3
Ref. B3324032	32 mm	3
Ref. B3324034	34 mm	3

Tornillo canulado para pie derecho
(Rosca a izquierda)
Cannulated screw for right foot
(Left-hand threaded)



	Longitud Length	Unidades Units
Ref. B3324110	10 mm	3
Ref. B3324112	12 mm	3
Ref. B3324114	14 mm	3
Ref. B3324116	16 mm	7
Ref. B3324118	18 mm	7
Ref. B3324120	20 mm	7
Ref. B3324122	22 mm	3
Ref. B3324124	24 mm	3
Ref. B3324126	26 mm	3
Ref. B3324128	28 mm	3
Ref. B3324130	30 mm	3
Ref. B3324132	32 mm	3
Ref. B3324134	34 mm	3

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IOSSUTQ01 / 07-2015 © Rev.1